DILABO, S.A.

CL Zacatena, 19
PG Industrial La Vega
13160 Torralba de Calatrava
Ciudad Real
Teléfono: 926 215 572
Fax: 926 215 612
Email: dilabo@dilabo.com
www.dilabo.com

MICROSCOPIOS METALGRÁFICOS
METALLURGICAL MICROSCOPES

Referencia/Code: 50170007

Referencia/Code: 50170008

Referencia / Code	50170007	50170008	50170072	50170082				
Modelo / Model	177	178	177/2	178/2				
Cabezal	Binocular	Triocular	Binocular	Triocular				
Head	Binocular	Trinocular	Binocular	Trinocular				
Oculares / Eyepieces	WF10x							
Revolcer / Nosepiece	Cuádruple							
Objetivos / Objectives	Acromáticos / Achromatic: 4x (A.N.: 0.10) 10x (A.N.: 0.25) 40x (R) (A.N.: 0.65) 60x (R) (A.N.: 0.85) Plano acromáticos: / Plan achromatic: 4x (A.N.: 0.10) 10x (A.N.: 0.25) 40x (R) (A.N.: 0.65) 60x (R) (A.N.: 0.85)							
Platina					Cuadrada 150x180 mm; orificio central 65 mm Ø Discos reductores de 12 y 20 mm Ø			
Stage					Square 150x180 mm; central hole 65 mm Ø Adapters 12 and 20 mm Ø			
Sistema iluminación					Köhler			
Illumination system					Köhler			
Condensador / Condenser	Abbe, regulable horizontalmente / Abbe, horizontally adjustable							
Diafragmas / Diaphragms	De campo y apertura / Field and aperture							
Filtros / Filters	Verde, azul y amarillo / Green, blue and yellow							
Lámpara	Halógena 6 V, 20 W con regulador de intensidad							
Lamp	Halogen 6 V, 20 W with intensity adjustment							

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We increase our offer of Zuzi metallurgical microscope for the study of the surface of any kind of opaque samples such as metals, alloys or minerals by the incidence of a light beam that reflects on the sample and reaches the objective. Therefore, these microscopes are of wide use and application in industry, laboratories of materials, geology, etc.

Ideal for large samples

The inverted design makes of these microscopes especially intended to observation of large samples impossible to observe with an up-right microscope due to their size. Furthermore, they are provided with a square stage (150x180 mm) with a central hole of 65 mm diameter. Each equipment is also supplied with 2 adapters with a central hole of 20 and 12 mm respectively to be used with specimens of several sizes.

High technical features

The head, binocular or trinocular (depending on the model), is 30° inclined and presents interpupillary distance adjustment (55-75 mm) as well as dioptric correction mechanism in one of the eyepieces. The trinocular microscopes are provided with a third tube with a sliding rod to direct light, that allows the use of photographic or video cameras to record images of the different observations made, increasing the features of the equipment.

The sample focusing is performed through coaxial coarse and fine focus adjustment knobs located at both sides of the microscope and provided with a friction control ring and a focusing lower limit blocking lever to avoid the slide to accidentally hit the objectives. Moreover, the fine adjustment knobs are graduated.

Excellent illumination of the sample

Zuzi metallurgical microscopes are equipped with a Köhler illumination system that allows a better control of the light beam that reaches the sample. This illumination system includes an Abbe condenser, field and aperture diaphragms and an halogen lamp with intensity adjustment. Microscopes are also supplied with filter holder and blue, green and yellow filters.

The lamp presents a centring mechanism through two knobs to get a correct adjustment and thus assuring an optimum illumination and visualization of the samples.